

# HEDIS PRACTICE SETUP

## NAVIGATION

Medical Records > Medical Records Setup > DM/Quality > HEDIS > HEDIS Practice Setup

NOTE: You must have security to access this area.

## ADD MEASURES

Measures are listed in order by ID.

All measures display in **red** until setup is complete. Each measure needs at least one CPT, Provider and Insurance.

## SELECT YOUR CODES

Your practice is responsible for setting up and maintaining each measure.

The codes you are billing must match the codes entered here; otherwise, you will not receive a pop-up.

## AUTO ANSWER

Additional Filters > Auto Answer Options

Check "Auto Answer Active" and select where you would like the Auto Answer to occur. Choose whether you'd like a pop-up if the results can't be determined and set your filters.

NOTE: This feature cannot be used with all measures (e.g. Falls Risk).

## ADDITIONAL FILTERS

You must make sure you have the correct labs and CPT codes selected to ensure the pop-up and/or Auto Answer will be triggered.

You must also make sure the status or "Date To Use" matches the lab result (e.g. RcvdCp "Received Complete").

## POP-UPS

When set up correctly, MEDENT will prompt you to choose an applicable CPT II or G-code for a qualifying patient in the following areas: Progress Note, eSuperbill and/or Charge Entry.