



















# MEDENT's

## Medication Selection Screen

The patient's name, account number, date of birth, age, sex and weight information will display along the top of this screen. Any drugs that have been discontinued by the manufacturer will highlight in a light gray in the medication list.

 <p>Click the chevron and gear to choose which icons will display under the chevron.</p>	 <p>Click the camera to view an image of the drug.</p>
 <p>OTC indicates it is an Over The Counter drug.</p>	 <p>G indicates the drug is generic.</p>
 <p>Click to view the dosing info. Note, this will not display if the patient is younger than the minimum dosing info in the MediSpan file.</p>	 <p>Click to view the drug ingredients. Note, this will only display if there is ingredient information.</p>
 <p>A green check mark indicates the medication is on formulary.</p>	 <p>Two green check marks indicate the medication is on formulary and the level preferred.</p>
 <p>NF indicates the medication is not on formulary.</p>	 <p>An NF with an X through it means the medication is not covered.</p>
 <p>NA indicates the formulary status is unknown.</p>	 <p>Alt Rx shows if therapeutic alternative meds exist. Click it to see the alternative meds and choose one to prescribe. When used with Surescripts, the drugs will show in order of formulary status level with the highest first.</p>
 <p>C shows the drug coverage information such as age, gender and quantity limits, medical necessity and step therapy. When a drug is selected that has additional coverage information available, the info will carry to the detail screen and display in red above the drug name.</p>	 <p>A dollar sign displays next to a drug when there is copay information available for the selected coverage. When clicked, a window will open showing what information is available for copay. Note, when a flat, percentage or tier copay is available for a drug, it will show as such.</p>
 <p><b>C</b> shows what classifications are assigned to that med. Select ALL CLASSIFICATIONS or a specific classification to view other meds assigned to that classification.</p>	 <p>Click <b>B</b> to only select the name of the drug. This lets the user select a "blank" drug. The <b>Oth Presc</b> box will automatically be checked when this is selected.</p>
 <p>Click this to retrieve coverage information for the patient. If there is at least one active coverage, a <b>Rx History</b> button will display, which can be used to retrieve the patient's medication history.</p>	 <p>This checks the formulary status of the patient's current medication list. When there are multiple coverages, the prescriber can enter the coverage window to select the coverage they wish to display.</p>