



# So you want to set up Windows Hello to sign in.

## What's next?

medent supports Windows Hello, a secure and time-saving biometric authentication option available on compatible Windows devices.

Instead of typing a username and password, users can log in with just a glance or a touch using facial recognition or fingerprint scanning.




## To Enable Windows Hello


- Open Settings on your Windows device.




- Click Accounts.



- Click Sign-in Options.




- Here, you can set up facial and/or fingerprint recognition and a PIN with Windows Hello.



## Within medent

- From the Main Menu, under Medical Records, click **Preferences > My Preferences > My Biometric (Windows Hello) Login**.



- Click Enable.
- Look at your webcam or swipe your finger.

**NOTE:** A practice administrator must also enable biometrics in **Password Rules**. Under Practice Management, click **Setup > Users/Depts/Security > Security System**. Then, select Password Rules and set **Allow Biometrics** to **Yes**.



- Once you have set up Windows Hello, when you open medent, you will be prompted to use the last biometric method you used. If you want to use a username and password, click Cancel.



Questions? Submit a support request at [medent.com](https://medent.com).